

TIGD5 Antibody

Catalog No: #42884

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	TIGD5 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification.
Applications	IHC
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total TIGD5 protein.
Immunogen Description	Fusion protein of human TIGD5
Target Name	TIGD5
Other Names	FLJ14926; MGC44883; TIGD5; tigger transposable element derived 5; Tigger transposable element-derived protein 5
Accession No.	Swiss-Prot#: Q53EQ6Gene ID: 84948
Uniprot	Q53EQ6
GeneID	84948;
Concentration	1.1mg/ml
Formulation	Rabbit IgG in pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol.
Storage	Store at -20°C


Application Details

Immunohistochemistry: 1:100-1:200

Images



Immunohistochemical analysis of paraffin-embedded Human Lung cancer tissue using #42884 at dilution 1/100.



Immunohistochemical analysis of paraffin-embedded Human Colorectal cancer tissue using #42884 at dilution 1/100.

Background

The protein encoded by this gene belongs to the tigger subfamily of the pogo superfamily of DNA-mediated transposons in humans. These proteins are related to DNA transposons found in fungi and nematodes, and more distantly to the Tc1 and mariner transposases. They are also very similar to the major mammalian centromere protein B. The exact function of this gene is not known.?

Note: This product is for in vitro research use only