

Recombinant Human DPPA2

Catalog No: #GP10374



Package Size: #GP10374-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human DPPA2
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 92-298 amino acids of human developmental pluripotency associated 2
Target Name	developmental pluripotency associated 2
Other Names	CT100; PESCRG1; ECAT15-2
Accession No.	Swissprot:Q7Z7J5Gene Accession:BC018070
Uniprot	Q7Z7J5
GenID	151871;
Storage	-20~-80°C, pH 7.6 PBS

Background

This protein binds to target gene promoters, including NKX2-5 and SYCE1, but not GATA4, and may be involved in the maintenance of the active epigenetic status of these genes. Expressed in embryonic stem cells. No expression is seen in 5 months embryo, mesenchymal stem cells, embryonic fibrocytes and adult tissues. It contains 1 SAP domain.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only