

# Recombinant Human DNM3

Catalog No: #GP10376



Package Size: #GP10376-1 100ug

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human DNM3
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 256-555 amino acids of human dynamin 3
Target Name	dynamin 3
Other Names	Dyna III
Accession No.	Swissprot:Q9UQ16Gene Accession:BC064546
Uniprot	Q9UQ16
GeneID	26052;
Storage	-20~-80°C, pH 7.6 PBS

## Background

This gene encodes a member of a family of guanosine triphosphate (GTP)-binding proteins that associate with microtubules and are involved in vesicular transport. The encoded protein functions in the development of megakaryocytes. Alternative splicing results in multiple transcript variants.

## References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only