

Recombinant Human ECE1

Catalog No: #GP10380



Package Size: #GP10380-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human ECE1
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to C terminal 300 amino acids of human endothelin converting enzyme 1
Target Name	endothelin converting enzyme 1
Other Names	ECE
Accession No.	Swissprot:P42892Gene Accession:BC117256
Uniprot	P42892
GeneID	1889;
Storage	-20~-80°C, pH 7.6 PBS

Background

The protein encoded by this gene is involved in proteolytic processing of endothelin precursors to biologically active peptides. Mutations in this gene are associated with Hirschsprung disease, cardiac defects and autonomic dysfunction. Alternatively spliced transcript variants encoding different isoforms have been noted for this gene.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only