

Recombinant Human ENPP5

Catalog No: #GP10390



Package Size: #GP10390-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human ENPP5
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 25-325 amino acids of human ectonucleotide pyrophosphatase/phosphodiesterase 5 (putative)
Target Name	ectonucleotide pyrophosphatase/phosphodiesterase 5 (putative)
Other Names	NPP-5
Accession No.	Swissprot:Q9UJA9Gene Accession:BC027615
Uniprot	Q9UJA9
GeneID	59084;
Storage	-20~-80°C, pH 7.6 PBS

Background

ENPP5 is a single-pass membrane protein that belongs to the nucleotide pyrophosphatase/phosphodiesterase family. It may play a role in neuronal cell communication, and it lacks nucleotide pyrophosphatase and lysopholipase D activity.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only