

Recombinant Human EDN2

Catalog No: #GP10405



Package Size: #GP10405-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human EDN2
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 25-178 amino acids of human endothelin 2
Target Name	endothelin 2
Other Names	ET2, PPET2
Accession No.	Swissprot:P20800Gene Accession:BC034393
Uniprot	P20800
GeneID	1907;
Storage	-20~-80°C, pH 7.6 PBS

Background

Endothelin 2, also known as ET-2, is a protein that in humans is encoded by the EDN2 gene. This gene encodes a member of the endothelin protein family of secretory vasoconstrictive peptides. The preproprotein is processed to a short mature form which functions as a ligand for the endothelin receptors that initiate intracellular signaling events. This gene product is involved in a wide range of biological processes, such as hypertension and ovulation.?

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only