

Recombinant Human FABP4

Catalog No: #GP10408



Package Size: #GP10408-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human FABP4
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 1-132 amino acids of human FABP4
Target Name	fatty acid binding protein 4
Other Names	aP2; ALBP; AFABP; A-FABP; HEL-S-104
Accession No.	Swissprot:P15090Gene Accession:BC003672
Uniprot	P15090
GeneID	2167;
Storage	-20~-80°C, pH 7.6 PBS

Background

FABP4 encodes the fatty acid binding protein found in adipocytes. Fatty acid binding proteins are a family of small, highly conserved, cytoplasmic proteins that bind long-chain fatty acids and other hydrophobic ligands. It is thought that FABPs roles include fatty acid uptake, transport, and metabolism.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only