

Recombinant Human FGL1

Catalog No: #GP10415



Package Size: #GP10415-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human FGL1
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 23-312 amino acids of human fibrinogen-like 1
Target Name	fibrinogen-like 1
Other Names	HFREP1; HP-041; LFIRE1; LFIRE-1
Accession No.	Swissprot:Q08830Gene Accession:BC007047
Uniprot	Q08830
GeneID	2267;
Storage	-20~-80°C, pH 7.6 PBS

Background

Fibrinogen-like 1 is a member of the fibrinogen family. This protein is homologous to the carboxy terminus of the fibrinogen beta- and gamma- subunits which contains the four conserved cysteines of fibrinogens and fibrinogen related proteins. However, this protein lacks the platelet-binding site, cross-linking region and a thrombin-sensitive site which are necessary for fibrin clot formation. This protein may play a role in the development of hepatocellular carcinomas. Four alternatively spliced transcript variants encoding the same protein exist for this gene.?

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only