

Recombinant Human FOXR2

Catalog No: #GP10418



Package Size: #GP10418-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human FOXR2
Brief Description	Recombinant Protein
Immunogen Description	Full length fusion protein
Target Name	forkhead box R2
Other Names	FOXN6
Accession No.	Swissprot:Q6PJQ5Gene Accession:BC012934
Uniprot	Q6PJQ5
GeneID	139628;
Storage	-20~-80°C, pH 7.6 PBS

Background

FOXR2 (Forkhead box protein R2), also known as FOXN6 (Forkhead box protein N6), is a 311 amino acid nuclear protein that belongs to the FOX family and contains one fork-head DNA-binding domain. Like other members of the FOX family, FOXR2 is not abundantly expressed in normal tissues but is found at very high levels in breast cancer and primary cancer cell lines and is thought to act as a transcription factor.?

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only