

Recombinant Human TBX20

Catalog No: #GP10420

Package Size: #GP10420-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human TBX20
Brief Description	Recombinant Protein
Immunogen Description	Full length fusion protein
Target Name	T-box 20
Other Names	ASD4
Accession No.	Swissprot:Q9UMR3Gene Accession:BC120946
Uniprot	Q9UMR3
GeneID	57057;
Storage	-20~-80°C, pH 7.6 PBS

Background

This gene encodes a T-box family member. The T-box family members share a common DNA binding domain, termed the T-box, and they are transcription factors involved in the regulation of developmental processes. This gene is essential for heart development. Mutations in this gene are associated with diverse cardiac pathologies, including defects in septation, valvulogenesis and cardiomyopathy. Alternatively spliced transcript variants encoding different isoforms have been found for this gene.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only