

Recombinant Human TBX6

Catalog No: #GP10421



Package Size: #GP10421-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human TBX6
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 274-436 amino acids of human TBX6
Target Name	T-box 6
Other Names	SCDO5
Accession No.	Swissprot:O95947Gene Accession:BC026031
Uniprot	O95947
GeneID	6911;
Storage	-20~-80°C, pH 7.6 PBS

Background

This gene is a member of a phylogenetically conserved family of genes that share a common DNA-binding domain, the T-box. T-box genes encode transcription factors involved in the regulation of developmental processes. Knockout studies in mice indicate that this gene is important for specification of paraxial mesoderm structures.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only