

## Recombinant Human GADD45G

Catalog No: #GP10430



Package Size: #GP10430-1 100ug

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human GADD45G
Brief Description	Recombinant Protein
Immunogen Description	Full length fusion protein
Target Name	growth arrest and DNA-damage-inducible, gamma
Other Names	CR6, DDIT2, GRP17, GADD45gamma
Accession No.	Swissprot:O95257Gene Accession:BC000465
Uniprot	O95257
GeneID	10912;
Storage	-20~-80°C, pH 7.6 PBS

## Background

This gene is a member of a group of genes whose transcript levels are increased following stressful growth arrest conditions and treatment with DNA-damaging agents. The protein encoded by this gene responds to environmental stresses by mediating activation of the p38/JNK pathway via MTK1/MEKK4 kinase. The GADD45G is highly expressed in placenta.

## References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only