

# Recombinant Human LGALS8

Catalog No: #GP10432



Package Size: #GP10432-1 100ug

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human LGALS8
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 187-317 amino acids of human lectin, galactoside-binding, soluble, 8
Target Name	lectin, galactoside-binding, soluble, 8
Other Names	Gal-8, PCTA1, PCTA-1, Po66-CBP
Accession No.	Swissprot:O00214Gene Accession:BC015818
Uniprot	O00214
GenID	3964;
Storage	-20~-80°C, pH 7.6 PBS

## Background

This gene encodes a member of the galectin family. Galectins are beta-galactoside-binding animal lectins with conserved carbohydrate recognition domains. The galectins have been implicated in many essential functions including development, differentiation, cell-cell adhesion, cell-matrix interaction, growth regulation, apoptosis, and RNA splicing. This gene is widely expressed in tumoral tissues and seems to be involved in integrin-like cell interactions. Alternatively spliced transcript variants encoding different isoforms have been identified.?

## References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only