

# Recombinant Human LGALS2

Catalog No: #GP10433



Package Size: #GP10433-1 100ug

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human LGALS2
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 1-132 amino acids of human LGALS2
Target Name	galectin 2
Other Names	HL14
Accession No.	Swissprot:P05162Gene Accession:BC059782
Uniprot	P05162
GeneID	3957;
Storage	-20~-80°C, pH 7.6 PBS

## Background

The protein encoded by this gene is a soluble beta-galactoside binding lectin. The encoded protein is found as a homodimer and can bind to lymphotoxin-alpha. A single nucleotide polymorphism in an intron of this gene can alter the transcriptional level of the protein, with a resultant increased risk of myocardial infarction.

## References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only