

Recombinant Human GAS2

Catalog No: #GP10435



Package Size: #GP10435-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human GAS2
Brief Description	Recombinant Protein
Immunogen Description	Full length fusion protein
Target Name	growth arrest-specific 2
Accession No.	Swissprot:O43903Gene Accession:BC013326
Uniprot	O43903
GeneID	2620;
Storage	-20~-80°C, pH 7.6 PBS

Background

The protein encoded by this gene is a caspase-3 substrate that plays a role in regulating microfilament and cell shape changes during apoptosis. It can also modulate cell susceptibility to p53-dependent apoptosis by inhibiting calpain activity. Multiple alternatively spliced variants, encoding the same protein, have been identified.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only