

# Recombinant Human CXCL3

Catalog No: #GP10437



Package Size: #GP10437-1 100ug

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human CXCL3
Brief Description	Recombinant Protein
Immunogen Description	Full length fusion protein
Target Name	chemokine (C-X-C motif) ligand 3
Other Names	GRO3; GROg; MIP2B; SCYB3; MIP-2b; CINC-2b
Accession No.	Swissprot:P19876Gene Accession:BC016308
Uniprot	P19876
GeneID	2921;
Storage	-20~-80°C, pH 7.6 PBS

## Background

Chemokine (C-X-C motif) ligand 3 (CXCL3) is a small cytokine belonging to the CXC chemokine family that is also known as GRO3 oncogene (GRO3), GRO protein gamma (GROg) and macrophage inflammatory protein-2-beta (MIP2b). CXCL3 controls migration and adhesion of monocytes and mediates its effects on its target cell by interacting with a cell surface chemokine receptor called CXCR2. More recently, it has been shown that Cxcl3 regulates cell autonomously the migration of the precursors of cerebellar granule neurons toward the internal layers of cerebellum, during the morphogenesis of cerebellum.

## References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only