

# Recombinant Human GDF3

Catalog No: #GP10441



Package Size: #GP10441-1 100ug

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human GDF3
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 251-364 amino acids of human growth differentiation factor 3
Target Name	growth differentiation factor 3
Other Names	KFS3, MCOP7, MCOPCB6
Accession No.	Swissprot:Q9NR23Gene Accession:BC030959
Uniprot	Q9NR23
GenelD	9573;
Storage	-20~-80°C, pH 7.6 PBS

## Background

The protein encoded by this gene is a member of the bone morphogenetic protein (BMP) family and the TGF-beta superfamily. This group of proteins is characterized by a polybasic proteolytic processing site which is cleaved to produce a mature protein containing seven conserved cysteine residues. The members of this family are regulators of cell growth and differentiation in both embryonic and adult tissues.

## References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only