

RABL2B Antibody

Catalog No: #42892

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	RABL2B Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification.
Applications	IHC
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total RABL2B protein.
Immunogen Type	protein
Immunogen Description	Fusion protein of human RABL2B
Target Name	RABL2B
Other Names	FLJ93981; FLJ98216; RAB, member of RAS oncogene family-like 2B; Rab-like protein 2B; RABL2B; RBL2B
Accession No.	Swiss-Prot#: Q9UNT1 Gene ID: 11158
Uniprot	Q9UNT1
GeneID	11158;
Concentration	1.5mg/ml
Formulation	Rabbit IgG in pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol.
Storage	Store at -20°C

Application Details

Immunohistochemistry: 1:25-1:100

Images



Immunohistochemical analysis of paraffin-embedded Human breast cancer tissue using #42892 at dilution 1/25,

Background

The RABL2B protein is a member of the RAB gene family which belongs to the RAS GTPase superfamily. RABL2B is located within a subtelomeric region of 22q13.3. Multiple alternatively spliced transcript variants encoding several different isoforms have been found for this gene.

Note: This product is for in vitro research use only