

Recombinant Human GPC6

Catalog No: #GP10450



Package Size: #GP10450-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human GPC6
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 229-529 amino acids of human glypican 6
Target Name	glypican 6
Other Names	OMIMD1
Accession No.	Swissprot:Q9Y625Gene Accession:BC106947
Uniprot	Q9Y625
GeneID	10082;
Storage	-20~-80°C, pH 7.6 PBS

Background

The glypicans comprise a family of glycosylphosphatidylinositol-anchored heparan sulfate proteoglycans, and they have been implicated in the control of cell growth and cell division. The glypican encoded by this gene is a putative cell surface coreceptor for growth factors, extracellular matrix proteins, proteases and anti-proteases.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only