

# Recombinant Human GRB7

Catalog No: #GP10452



Package Size: #GP10452-1 100ug

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human GRB7
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 232-532 amino acids of human growth factor receptor-bound protein 7
Target Name	growth factor receptor-bound protein 7
Accession No.	Swissprot:Q14451Gene Accession:BC006535
Uniprot	Q14451
GeneID	2886;
Storage	-20~-80°C, pH 7.6 PBS

## Background

The product of this gene belongs to a small family of adapter proteins that are known to interact with a number of receptor tyrosine kinases and signaling molecules. This gene encodes a growth factor receptor-binding protein that interacts with epidermal growth factor receptor (EGFR) and ephrin receptors. The protein plays a role in the integrin signaling pathway and cell migration by binding with focal adhesion kinase (FAK). Several transcript variants encoding two different isoforms have been found for this gene.

## References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only