

# Recombinant Human SPINT1

Catalog No: #GP10464



Package Size: #GP10464-1 100ug

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human SPINT1
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 36-140 amino acids of human serine peptidase inhibitor, Kunitz type 1
Target Name	serine peptidase inhibitor, Kunitz type 1
Other Names	HAI, HAI1, MANSC2
Accession No.	Swissprot:O43278Gene Accession:BC004140
Uniprot	O43278
GenID	6692;
Storage	-20~-80°C, pH 7.6 PBS

## Background

The protein encoded by this gene is a member of the Kunitz family of serine protease inhibitors. The protein is a potent inhibitor specific for HGF activator and is thought to be involved in the regulation of the proteolytic activation of HGF in injured tissues. Alternative splicing results in multiple variants encoding different isoforms.

## References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only