

Recombinant Human SCARB1

Catalog No: #GP10467



Package Size: #GP10467-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human SCARB1
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 465-552 amino acids of human scavenger receptor class B, member 1
Target Name	scavenger receptor class B, member 1
Other Names	CLA1; SRB1; CLA-1; SR-BI; CD36L1; HDLQTL6
Accession No.	Swissprot:Q8WTV0Gene Accession:BC022087
Uniprot	Q8WTV0
GenID	949;
Storage	-20~-80°C, pH 7.6 PBS

Background

The protein encoded by this gene is a plasma membrane receptor for high density lipoprotein cholesterol (HDL). The encoded protein mediates cholesterol transfer to and from HDL. In addition, this protein is a receptor for hepatitis C virus glycoprotein E2. Two transcript variants encoding different isoforms have been found for this gene.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only