

Recombinant Human HOXA6

Catalog No: #GP10480



Package Size: #GP10480-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human HOXA6
Brief Description	Recombinant Protein
Immunogen Description	Full length fusion protein
Target Name	homeobox A6
Other Names	HOX1, HOX1B, HOX1.2
Accession No.	Swissprot:P31267Gene Accession:BC069497
Uniprot	P31267
GeneID	3203;
Storage	-20~-80°C, pH 7.6 PBS

Background

In vertebrates, the genes encoding the class of transcription factors called homeobox genes are found in clusters named A, B, C, and D on four separate chromosomes. Expression of these proteins is spatially and temporally regulated during embryonic development. This gene is part of the A cluster on chromosome 7 and encodes a DNA-binding transcription factor which may regulate gene expression, morphogenesis, and differentiation.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only