

Recombinant Human HSF2

Catalog No: #GP10486



Package Size: #GP10486-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human HSF2
Brief Description	Recombinant Protein
Immunogen Description	Full length fusion protein
Target Name	heat shock transcription factor 2
Other Names	HSF 2, HSTF 2
Accession No.	Swissprot:Q03933Gene Accession:BC005329
Uniprot	Q03933
GeneID	3298;
Storage	-20~-80°C, pH 7.6 PBS

Background

The protein encoded by this gene belongs to the HSF family of transcription factors that bind specifically to the heat-shock promoter element and activate transcription. Heat shock transcription factors activate heat-shock response genes under conditions of heat or other stresses. Alternatively spliced transcript variants encoding different isoforms have been found for this gene.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only