

# Recombinant Human IFNG

Catalog No: #GP10499



Package Size: #GP10499-1 100ug

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human IFNG
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 24-161 amino acids of human IFNG
Target Name	interferon, gamma
Other Names	IFG; IFI
Accession No.	Swissprot:P01579Gene Accession:BC070256
Uniprot	P01579
GeneID	3458;
Storage	-20~-80°C, pH 7.6 PBS

## Background

This gene encodes a member of the type II interferon family. The protein encoded is a soluble cytokine with antiviral, immunoregulatory and anti-tumor properties and is a potent activator of macrophages. Mutations in this gene are associated with aplastic anemia.

## References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only