

# Recombinant Human IGFBP4

Catalog No: #GP10500



Package Size: #GP10500-1 100ug

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human IGFBP4
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 22-258 amino acids of human insulin-like growth factor binding protein 4
Target Name	insulin-like growth factor binding protein 4
Other Names	BP-4; IBP4; IGFBP-4; HT29-IGFBP
Accession No.	Swissprot:P22692Gene Accession:BC016041
Uniprot	P22692
GenID	3487;
Storage	-20~-80°C, pH 7.6 PBS

## Background

This gene is a member of the insulin-like growth factor binding protein (IGFBP) family and encodes a protein with an IGFBP domain and a thyroglobulin type-I domain. The protein binds both insulin-like growth factors (IGFs) I and II and circulates in the plasma in both glycosylated and non-glycosylated forms. Binding of this protein prolongs the half-life of the IGFs and alters their interaction with cell surface receptors.

## References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only