

Recombinant Human IGFBP5

Catalog No: #GP10501



Package Size: #GP10501-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human IGFBP5
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 21-272 amino acids of human insulin-like growth factor binding protein 5
Target Name	insulin-like growth factor binding protein 5
Other Names	IBP5
Accession No.	Swissprot:P24593Gene Accession:BC011453
Uniprot	P24593
GenID	3488;
Storage	-20~-80°C, pH 7.6 PBS

Background

Insulin-like growth factor-binding protein 5 is a protein that in humans is encoded by the IGFBP5 gene. IGF-binding proteins prolong the half-life of the IGFs and have been shown to either inhibit or stimulate the growth promoting effects of the IGFs on cell culture. They alter the interaction of IGFs with their cell surface receptors.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only