

# Recombinant Human IL3

Catalog No: #GP10511



Package Size: #GP10511-1 100ug

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human IL3
Brief Description	Recombinant Protein
Immunogen Description	Full length fusion protein
Target Name	Interleukin 3 (colony-stimulating factor, multiple)
Other Names	IL-3; MCGF; MULTI-CSF
Accession No.	Swissprot:P08700Gene Accession:BC069472
Uniprot	P08700
GeneID	3562;
Storage	-20~-80°C, pH 7.6 PBS

## Background

The protein encoded by this gene is a potent growth promoting cytokine. This cytokine is capable of supporting the proliferation of a broad range of hematopoietic cell types. It is involved in a variety of cell activities such as cell growth, differentiation and apoptosis. This cytokine has been shown to also possess neurotrophic activity, and it may be associated with neurologic disorders.

## References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only