

Recombinant Human IL18

Catalog No: #GP10514



Package Size: #GP10514-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human IL18
Brief Description	Recombinant Protein
Immunogen Description	Full length fusion protein
Target Name	interleukin 18
Other Names	IGIF, IL-18, IL-1g, IL1F4
Accession No.	Swissprot:Q14116Gene Accession:BC007007
Uniprot	Q14116
GeneID	3606;
Storage	-20~-80°C, pH 7.6 PBS

Background

The protein encoded by this gene is a proinflammatory cytokine that augments natural killer cell activity in spleen cells, and stimulates interferon gamma production in T-helper type I cells. Alternatively spliced transcript variants encoding different isoforms have been found for this gene

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only