

Recombinant Human IL20

Catalog No: #GP10515



Package Size: #GP10515-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human IL20
Brief Description	Recombinant Protein
Immunogen Description	Full length fusion protein
Target Name	interleukin 20
Other Names	IL-20; IL10D; ZCYTO10
Accession No.	Swissprot:Q9NYY1Gene Accession:BC069311
Uniprot	Q9NYY1
GeneID	50604;
Storage	-20~-80°C, pH 7.6 PBS

Background

The protein encoded by this gene is a cytokine structurally related to interleukin 10 (IL10). This cytokine has been shown to transduce its signal through signal transducer and activator of transcription 3 (STAT3) in keratinocytes. A specific receptor for this cytokine is found to be expressed in skin and upregulated dramatically in psoriatic skin, suggesting a role for this protein in epidermal function and psoriasis.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only