

## Recombinant Human ING2

Catalog No: #GP10519



Package Size: #GP10519-1 100ug

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human ING2
Brief Description	Recombinant Protein
Immunogen Description	Full length fusion protein
Target Name	inhibitor of growth family, member 2
Other Names	ING1L; p33ING2
Accession No.	Swissprot:Q9H160Gene Accession:BC030128
Uniprot	Q9H160
GeneID	3622;
Storage	-20~-80°C, pH 7.6 PBS

## Background

This gene is a member of the inhibitor of growth (ING) family. Members of the ING family associate with and modulate the activity of histone acetyltransferase (HAT) and histone deacetylase (HDAC) complexes and function in DNA repair and apoptosis

## References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only