

# Recombinant Human INHA

Catalog No: #GP10520



Package Size: #GP10520-1 100ug

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human INHA
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 233-366 amino acids of human inhibin, alpha
Target Name	inhibin, alpha
Accession No.	Swissprot:P05111Gene Accession:BC006391
Uniprot	P05111
GeneID	3623;
Storage	-20~-80°C, pH 7.6 PBS

## Background

This gene encodes the alpha subunit of inhibins A and B protein complexes. These complexes negatively regulate follicle stimulating hormone secretion from the pituitary gland. Inhibins have also been implicated in regulating numerous cellular processes including cell proliferation, apoptosis, immune response and hormone secretion. [

## References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only