

Recombinant Human ITGB3

Catalog No: #GP10523



Package Size: #GP10523-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|-----------------------|--|
| Product Name | Recombinant Human ITGB3 |
| Brief Description | Recombinant Protein |
| Immunogen Description | Fusion protein corresponding to a region derived from 418-718 amino acids of human integrin, beta 3 (platelet glycoprotein IIIa, antigen CD61) |
| Target Name | integrin, beta 3 (platelet glycoprotein IIIa, antigen CD61) |
| Other Names | GT; CD61; GP3A; BDPLT2; GPIIIa; BDPLT16 |
| Accession No. | Swissprot:P05106Gene Accession:BC127666 |
| Uniprot | P05106 |
| GenelD | 3690; |
| Storage | -20~-80°C, pH 7.6 PBS |

Background

The ITGB3 protein product is the integrin beta chain beta 3. Integrins are integral cell-surface proteins composed of an alpha chain and a beta chain. A given chain may combine with multiple partners resulting in different integrins. Integrin beta 3 is found along with the alpha IIb chain in platelets. Integrins are known to participate in cell adhesion as well as cell-surface mediated signalling.?

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only