

Recombinant Human IRAK3

Catalog No: #GP10535



Package Size: #GP10535-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human IRAK3
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 165-452 amino acids of human IRAK3
Target Name	interleukin-1 receptor-associated kinase 3
Other Names	ASRT5; IRAKM
Accession No.	Swissprot:Q9Y616Gene Accession:BC057800
Uniprot	Q9Y616
GeneID	11213;
Storage	-20~-80°C, pH 7.6 PBS

Background

This gene encodes a member of the interleukin-1 receptor-associated kinase protein family. Members of this family are essential components of the Toll/IL-R immune signal transduction pathways. This protein is primarily expressed in monocytes and macrophages and functions as a negative regulator of Toll-like receptor signaling. Mutations in this gene are associated with a susceptibility to asthma. Alternate splicing results in multiple transcript variants.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only