

Recombinant Human KLF11

Catalog No: #GP10546



Package Size: #GP10546-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human KLF11
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to C terminal 300 amino acids of human Kruppel-like factor 11
Target Name	Kruppel-like factor 11
Other Names	FKLF; FKLF1; MODY7; TIEG2; Tieg3
Accession No.	Swissprot:O14901Gene Accession:BC063286
Uniprot	O14901
GeneID	8462;
Storage	-20~-80°C, pH 7.6 PBS

Background

The protein encoded by this gene is a zinc finger transcription factor that binds to SP1-like sequences in epsilon- and gamma-globin gene promoters. This binding inhibits cell growth and causes apoptosis. Defects in this gene are a cause of maturity-onset diabetes of the young type 7 (MODY7). Three transcript variants encoding two different isoforms have been found for this gene.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only