

Recombinant Human KLF15

Catalog No: #GP10547



Package Size: #GP10547-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human KLF15
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to N terminal 300 amino acids of human Kruppel-like factor 15
Target Name	Kruppel-like factor 15
Other Names	KKLF
Accession No.	Swissprot:Q9UIH9Gene Accession:BC036733
Uniprot	Q9UIH9
GeneID	28999;
Storage	-20~-80°C, pH 7.6 PBS

Background

Kruppel-like factor 15 is a protein that in humans is encoded by the KLF15 gene in the Kruppel-like factor family. KLF15, a kidney-enriched Kruppel-like factor, is a transcriptional activator that binds the CLCNKA promoter. KLF15, which is a nuclear protein, is expressed primarily in liver, heart, skeletal muscle and kidney tissues but is not detected in lymphoid tissues or bone marrow. It is an important regulator of GLUT4 in both adipose and muscle tissues.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only