

Recombinant Human KLK1

Catalog No: #GP10549



Package Size: #GP10549-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human KLK1
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 19-262 amino acids of human kallikrein 1
Target Name	kallikrein 1
Other Names	hK1; KLKR; Klk6
Accession No.	Swissprot:P06870Gene Accession:BC005313
Uniprot	P06870
GeneID	3816;
Storage	-20~-80°C, pH 7.6 PBS

Background

Kallikreins are a subgroup of serine proteases having diverse physiological functions. Growing evidence suggests that many kallikreins are implicated in carcinogenesis and some have potential as novel cancer and other disease biomarkers. This gene is one of the fifteen kallikrein subfamily members located in a cluster on chromosome 19. This protein is functionally conserved in its capacity to release the vasoactive peptide, Lys-bradykinin, from low molecular weight kininogen.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only