

Recombinant Human LCAT

Catalog No: #GP10557



Package Size: #GP10557-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human LCAT
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 141-440 amino acids of human LCAT
Target Name	lecithin-cholesterol acyltransferase
Accession No.	Swissprot:P04180Gene Accession:BC014781
Uniprot	P04180
GeneID	3931;
Storage	-20~-80°C, pH 7.6 PBS

Background

This gene encodes the extracellular cholesterol esterifying enzyme, lecithin-cholesterol acyltransferase. The esterification of cholesterol is required for cholesterol transport. Mutations in this gene have been found to cause fish-eye disease as well as LCAT deficiency.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only