

Recombinant Human LGR5

Catalog No: #GP10559



Package Size: #GP10559-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human LGR5
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 22-320 amino acids of human leucine-rich repeat containing G protein-coupled receptor 5
Target Name	leucine-rich repeat containing G protein-coupled receptor 5
Other Names	FEX, HG38, GPR49, GPR67, GRP49
Accession No.	Swissprot:O75473Gene Accession:BC096324
Uniprot	O75473
GenelD	8549;
Storage	-20~-80°C, pH 7.6 PBS

Background

This gene encodes Orphan receptor. Stem cell marker of the intestinal epithelium and the hair follicle. Target gene of Wnt signaling.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only