

Recombinant Human LRP12

Catalog No: #GP10564



Package Size: #GP10564-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human LRP12
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 514-859 amino acids of human low density lipoprotein receptor-related protein 12
Target Name	low density lipoprotein receptor-related protein 12
Other Names	ST7
Accession No.	Swissprot:Q9Y561Gene Accession:BC032109
Uniprot	Q9Y561
GeneID	29967;
Storage	-20~-80°C, pH 7.6 PBS

Background

This gene encodes a member of the low-density lipoprotein receptor related protein family. The product of this gene is a transmembrane protein that is differentially expressed in many cancer cells. Alternate splicing results in multiple transcript variants.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only