

Recombinant Human LYVE1

Catalog No: #GP10568



Package Size: #GP10568-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human LYVE1
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 20-238 amino acids of human LYVE1
Target Name	lymphatic vessel endothelial hyaluronan receptor 1
Other Names	HAR, XLKD1, LYVE-1, CRSBP-1
Accession No.	Swissprot:Q9Y5Y7 Gene Accession:BC026231
Uniprot	Q9Y5Y7
GeneID	10894;
Storage	-20~-80°C, pH 7.6 PBS

Background

This gene encodes a type I integral membrane glycoprotein. The encoded protein acts as a receptor and binds to both soluble and immobilized hyaluronan. This protein may function in lymphatic hyaluronan transport and have a role in tumor metastasis.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only