

Recombinant Human MATN3

Catalog No: #GP10578



Package Size: #GP10578-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human MATN3
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to C terminal 300 amino acids of human matrilin 3
Target Name	matrilin 3
Other Names	HOA; OS2; EDM5; DIPOA; OADIP
Accession No.	Swissprot:O15232Gene Accession:BC139907
Uniprot	O15232
GeneID	4148;
Storage	-20~-80°C, pH 7.6 PBS

Background

This gene encodes a member of von Willebrand factor A domain containing protein family. This family of proteins is thought to be involved in the formation of filamentous networks in the extracellular matrices of various tissues. This protein contains two von Willebrand factor A domains; it is present in the cartilage extracellular matrix and has a role in the development and homeostasis of cartilage and bone. Mutations in this gene result in multiple epiphyseal dysplasia.?

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only