

# Recombinant Human TMPRSS7

Catalog No: #GP10581



Package Size: #GP10581-1 100ug

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human TMPRSS7
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to C terminal 300 amino acids of human transmembrane protease, serine 7
Target Name	transmembrane protease, serine 7
Accession No.	Swissprot:Q7RTY8Gene Accession:BC117322
Uniprot	Q7RTY8
GeneID	344805;
Storage	-20~-80°C, pH 7.6 PBS

## Background

Serine protease which preferentially hydrolyzes peptides with Arg at the P1 position. Expressed in brain, ovary, testis, salivary gland, trachea and lung.

## References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only