

Recombinant Human MMP28

Catalog No: #GP10602



Package Size: #GP10602-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human MMP28
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to C terminal 300 amino acids of human matrix metalloproteinase 28
Target Name	matrix metalloproteinase 28
Other Names	MM28, MMP25, MMP-28, EPILYSIN
Accession No.	Swissprot:Q9H239Gene Accession:BC002631
Uniprot	Q9H239
GeneID	79148;
Storage	-20~-80°C, pH 7.6 PBS

Background

Proteins of the matrix metalloproteinase (MMP) family are involved in the breakdown of extracellular matrix for both normal physiological processes, such as embryonic development, reproduction and tissue remodeling, and disease processes, such as asthma and metastasis. This gene encodes a secreted enzyme that degrades casein. Its expression pattern suggests that it plays a role in tissue homeostasis and in wound repair. Transcript variants encoding different isoforms have been described.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only