

Recombinant Human MYOZ1

Catalog No: #GP10614



Package Size: #GP10614-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human MYOZ1
Brief Description	Recombinant Protein
Immunogen Description	Full length fusion protein
Target Name	myozenin 1
Other Names	CS-2, FATZ, MYOZ
Accession No.	Swissprot:Q9NP98Gene Accession:BC025753
Uniprot	Q9NP98
GeneID	58529;
Storage	-20~-80°C, pH 7.6 PBS

Background

The protein encoded by this gene is primarily expressed in the skeletal muscle, and belongs to the myozenin family. Members of this family function as calcineurin-interacting proteins that help tether calcineurin to the sarcomere of cardiac and skeletal muscle. They play an important role in modulation of calcineurin signaling.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only