

## Recombinant Human MYL6

Catalog No: #GP10619



Package Size: #GP10619-1 100ug

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human MYL6
Brief Description	Recombinant Protein
Immunogen Description	Full length fusion protein
Target Name	myosin, light chain 6, alkali, smooth muscle and non-muscle
Other Names	LC17; ESMLC; LC17A; LC17B; MLC-3; MLC1SM; MLC3NM; MLC3SM; LC17-GI; LC17-NM
Accession No.	Swissprot:P60660Gene Accession:BC017455
Uniprot	P60660
GeneID	4637;
Storage	-20~-80°C, pH 7.6 PBS

## Background

Myosin is a hexameric ATPase cellular motor protein. It is composed of two heavy chains, two nonphosphorylatable alkali light chains, and two phosphorylatable regulatory light chains. This gene encodes a myosin alkali light chain that is expressed in smooth muscle and non-muscle tissues. Genomic sequences representing several pseudogenes have been described and two transcript variants encoding different isoforms have been identified for this gene.

## References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only