

Recombinant Human NAP1L1

Catalog No: #GP10622



Package Size: #GP10622-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human NAP1L1
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 1-300 amino acids of human NAP1L1
Target Name	nucleosome assembly protein 1-like 1
Other Names	NRP, NAP1, NAP1L
Accession No.	Swissprot:P55209Gene Accession:BC002387
Uniprot	P55209
GeneID	4673;
Storage	-20~-80°C, pH 7.6 PBS

Background

This gene encodes a member of the nucleosome assembly protein (NAP) family. This protein participates in DNA replication and may play a role in modulating chromatin formation and contribute to the regulation of cell proliferation. Alternative splicing of this gene results in several transcript variants; however, not all have been fully described.?

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only