

# Recombinant Human NGB

Catalog No: #GP10625



Package Size: #GP10625-1 100ug

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human NGB
Brief Description	Recombinant Protein
Immunogen Description	Full length fusion protein
Target Name	neuroglobin
Accession No.	Swissprot:Q9NPG2Gene Accession:BC032509
Uniprot	Q9NPG2
GeneID	58157;
Storage	-20~-80°C, pH 7.6 PBS

## Background

This gene encodes an oxygen-binding protein that is distantly related to members of the globin gene family. It is highly conserved among other vertebrates. It is expressed in the central and peripheral nervous system where it may be involved in increasing oxygen availability and providing protection under hypoxic/ischemic conditions.?

## References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only