

Recombinant Human NEUROG3

Catalog No: #GP10631



Package Size: #GP10631-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human NEUROG3
Brief Description	Recombinant Protein
Immunogen Description	Full length fusion protein
Target Name	neurogenin 3
Other Names	ngn3; Atoh5; NGN-3; Math4B; bHLHa7
Accession No.	Swissprot:Q9Y4Z2Gene Accession:BC117488
Uniprot	Q9Y4Z2
GeneID	50674;
Storage	-20~-80°C, pH 7.6 PBS

Background

The protein encoded by this gene is a basic helix-loop-helix (bHLH) transcription factor involved in neurogenesis. The encoded protein likely acts as a heterodimer with another bHLH protein. Defects in this gene are a cause of congenital malabsorptive diarrhea 4 (DIAR4).

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only