

Recombinant Human NME5

Catalog No: #GP10640



Package Size: #GP10640-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human NME5
Brief Description	Recombinant Protein
Immunogen Description	Full length fusion protein
Target Name	NME/NM23 family member 5
Other Names	NM23H5; RSPH23; NM23-H5
Accession No.	Swissprot:P56597Gene Accession:BC026182
Uniprot	P56597
GeneID	8382;
Storage	-20~-80°C, pH 7.6 PBS

Background

The NME5 gene encodes a 212-amino acid polypeptide with homology to the NDP kinases. Northern blot analysis detected abundant expression of a 1.1-kb mRNA in testis. In situ hybridization to testis revealed NME5 mRNA in the seminiferous tubules in a pattern consistent with stage-specific expression in germinal cells during spermatogenesis.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only